

Benzoic Acid Tablets (Standard for Calorific Analysis) and Fuse Wire


Part No.	Description	Tablet	Unit of Sale
CA-208C	Benzoic Acid Tablets 1.0 gm (100/btl)	1.0 gm	100/btl
CA-208V	Benzoic Acid Tablets 1.0 gm (500/btl) (now sold as 5 x CA-208C at discounted price)	1.0 gm	500/btl
CA-0012C	Fuse Wire Roll (100 gm)		100 gms
CA-0012M	Fuse Wire Roll (1000 gm)		1000 gms
AR-2134C	Benzoic Acid Tablets 1.45 gm (100/btl)	1.45 gm	100/btl
AR-3403	Benzoic Acid POWDER (454 gms) for SPIKING		454 gms
AR-3414V	Benzoic Acid Tablets (500 x 0.2 gm) for SPIKING	0.20 gm	500/btl
ARA-4204	S.S. Crucibles for Calorimeters		10/pkt
ARA-8671	In-Line Filter for AC-350		

Sundry Reagents etc.


Part No.	Description	Unit of Sale
AR-060	Manganese Dioxide (coarse), absorbs SO ₂ (30 cc = 60 gm)	30 cc
AR-062	Manganese Dioxide (coarse), absorbs SO ₂ (250 cc = 500 gm)	250 cc
AR-081	Glass Wool	1 lb = 454 gms
AR-1608	Quartz Wool Wadding (coarse)	200 gms
AR-168	AlphaSolve (absorbs CO ₂)	250 gms
AR-177-50	Quartz Wool (fine)	50 gms
AR-611	Schutze's Reagent (CO to CO ₂)	100 gms
AR-636	Vanadium Pentoxide (O ₂ provider)	250 gms
CA-171	Anhydrone (absorbs H ₂ O, MgClO ₃) Discounted for purchase of large quantities	1 lb = 454 gms
CA-8379	Quartz Wool Strips	10/pkt
CA-CERAM/500	Autostatic Ceramic Adhesive (500gms)	500gms
CU2P207-H2O	Copper (II) Pyrophosphate Hydrate	500 gms
SILVGAUZE	Silver Gauze (absorbs SO ₂)	each

Tiles and Items for Ash Fusion Determination

Part No.	Description	Unit of Sale
92-2833-01	Gold Wire 0.5 mm diam 99.99% 305mm M.P. 1064 °C	305 mm / pk
92-2837-01	Nickel Wire 0.5mm diam 99.99% 305mm M.P. 1453 °C	305 mm / pk
AFT-42424	Ash Fusion Tile - 24 x 24 x 3 mm with 2 Bevelled Edges	100 / pk
AFT-62525	Ash Fusion Tile - Flat 25 x 25 x 6 mm	200 / pk
Pd-Wire	Palladium Wire 0.5 mm diam 99.95% 500mm M.P. 1552 °C	500 mm / pk



AFT-62525



AFT-42424

60015 **Marking Ink** for Ceramic Items - 15 ml.
For permanently marking laboratory porcelain and other ceramic ware. May be applied with pen, brush or toothpick. Marked surface must be heat-treated to become permanent. Use instructions on bottle. In widemouth bottle with screw cap.



010-070-020 **Marking Pen** - High Temperature Marker (to 1100°C)



Please see other pages for extensive range of Standards Australia ASCRM and Choice Analytical Coal and Coke Standards.



CHOICE Analytical Pty Ltd

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.au

Web : www.choiceanalytical.com.au