


## ICP Mass Spectrometer Skimmers and Sampler Cones

Choice Analytical is now offering high quality ICP mass spectrometer cones and skimmers. These items are supplied at realistic prices and are guaranteed to perform and last as good or better than original equipment manufacturers' (OEM) supplied parts. The cones and skimmers are CNC machined with high quality alloys and acute precision for quality performance. Consult the chart below for material selection and OEM part number cross reference.

Part No. OEM	In sequence of Choice Part No.	Description
--------------	-----------------------------------	-------------


### VG PLASMAQUAD / PLASMATRACE

3004661/4001690	APT-9001	Nickel Sampler
3004866/4001689	APT-9002	Nickel Mini-Skimmer
3200860/4200846	APT-9003	Nickel Micro-Skimmer (EI)*
3200262	APT-9004	Platinum Sampler 18 mm
—	APT-9005	Platinum Mini-Skimmer
3201101	APT-9006	Platinum Micro-Skimmer (EI)*
3202005	APT-9007	Platinum Plasmatrace Skimmer
3200456	APT-9008	Nickel Plasmatrace Skimmer
3201941	APT-9007B	Platinum Plasmatrace II Micro Skimmer
3201457	APT-9008B	Nickel Plasmatrace II Micro-Skimmer




### PERKIN-ELMER

N810-2101	APT-9009	Nickel Sampler 250/500
N810-2111	APT-9010	Nickel Skimmer 250/500
WE01-3087	APT-9011	Nickel Sampler 5000
N812-2000	APT-9012	Nickel Skimmer 5000
N810-2131	APT-9013	Platinum Sampler 250/500
N810-2141	APT-9014	Platinum Skimmer 250/500
WE01-3088	APT-9015	Platinum Sampler 5000
N812-2035	APT-9016	Platinum Skimmer 5000
WE02-1140	APT-9017	Nickel Sampler 6000
WE02-1137	APT-9018	Nickel Skimmer 6000
WE01-4815	APT-9019	Platinum Sampler 6000
WE01-4816	APT-9020	Platinum Skimmer 6000



### FINNIGAN SOLA

TS067014	APT-9021	Nickel Sampler
FS067021	APT-9022	Nickel Skimmer MKIII Enhanced
FS067019	APT-9023	Nickel Accelerator
TS067004	APT-9024	Platinum Sampler 18 mm insert
TS067016	APT-9025	Platinum Skimmer 7.2 mm insert



### AGILENT / HEWLETT PACKARD

G1820-65238	APT-9026	Nickel Sampler 25 mm
G1820-65265	APT-9027	Nickel Skimmer
G1820-65239	APT-9028	Platinum Sampler 10 mm
—	APT-9029	Platinum Sampler 18 mm
G1820-65237	APT-9030	Platinum Skimmer 7.2 mm



### THERMO JARREL ASH

13307600	APT-9031	Nickel Sampler
13307500	APT-9032	Nickel Skimmer
13555800	APT-9033	Platinum Sampler 18 mm
13582700	APT-9034	Platinum Skimmer 7.2 mm



### VARIAN

99-100698	APT-9035	Nickel Sampler (Cu Base with 1.5" Ni Insert)
99-100699	APT-9036	Nickel Skimmer
99-100700	APT-9037	Platinum Sampler 18 mm
99-100701	APT-9038	Platinum Skimmer 7.2 mm

\* - (EI) Enhanced Interface)

