

### Peristaltic Pump Tubing:-

- Tygon / PVC
- Coated PVC
- Viton
- Silicone

Suitable for use on:-

**AGILENT / HP**

ICP-MS



All tubing features **two or three colour-coded stops** for easy determination of size and is available in • **Standard Tygon / PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Silicone** • **Santoprene\* (creamy opaque)**. Tubing overall length is 16", with a 6" centre pumping section separated by two colour-coded stops, three stops (72 mm between stops) Our Viton Tubing (black opaque) has an overall length of 7" and also includes two or three stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID inches/mm	Agilent Part No.	Standard Tygon/PVC Tubing	Coated PVC Solvent Resistant Tubing	Santoprene Tubing *	Pack Size 12/pk	Silicone Tubing	Pack Size 6/pk
<b>TWO-STOP TYPE</b>								
White/White	0.040/1.02	5064-8028	AR-4311	-	-	12/pk	-	-
Orange/Blue	0.010/0.25	5064-8015	AR-4302	-	-	12/pk	-	-
<b>THREE-STOP TYPE</b>								
used on ICP-MS, FIAS and FIMS systems								
White/White/White	0.040/1.02	G1820-65215	AR-34311	-	-	12/pk	-	-
White/White/White	0.040/1.02	G1820-65217	-	-	-	-	AR-34311S	6/pk
Yellow/Blue/Yellow	0.060/1.52	G1820-65216	-	-	AR-34317SA	12/pk	-	-
Orange/Red/Orange	0.008/0.19	G1820-65220	AR-34301	-	-	12/pk	-	-
Orange/Blue/Orange	0.010/0.25	5064-8034	AR-34302	-	-	12/pk	-	-
Orange/Orange/Orange	0.035/0.89	5064-8014	AR-34309	-	-	12/pk	-	-

Many other sizes and types of Tubing and Pump Tubes are also available e.g. Viton.

Please enquire with full details of requirement.

**[\*] Santoprene = PharMed = Ismaprene**

**All are same material - creamy opaque**