


Listing of Equivalent Items available from **CHOICE ANALYTICAL** for use on PERKIN ELMER instrumentation and accessories



Sample Cups/Tubes used on Auto Samplers

In sequence of:

PERKIN-ELMER Part No.	Part No. CHOICE	Description	Unit of Sale
0047-1441	4092-314	15 mL Conical Bottom Tube, Polystyrene for PE AS-50/51	pkt/125
B008-7056	AAA-2001	2ml Polyethylene Conical Cups for PE AS-1, AS-40/60/70	pkt/1,000
B008-7056	AAA-2010	2ml Polyethylene Conical Cups for PE AS-1, AS-40/60/70	pkt/10,000
B011-9079	AAA-1001	2ml Polystyrene Conical Cups for PE AS-1, AS-40/60/70	pkt/1,000
B011-9079	AAA-1010	2ml Polystyrene Conical Cups for PE AS-1, AS-40/60/70	pkt/10,000
B012-9303	AAA-4001	4ml Polystyrene Conical Cups for PE AS-1, AS-60/70	pkt/1,000
B019-3233	4092-316	Centrifuge Tube, 15ml, Polypropylene, metal free for PE AS-90/91	pkt/500
—	4092-320	16-17mm x 100mm, Polypropylene Tube for PE AS-90	pkt/1,000
—	4092-321	16-17mm x 100mm, Polypropylene Tube for PE AS-90	pkt/10,000
—	4092-324	13 x 100mm Polystyrene Tube for PE	pkt/1,000
—	4092-327	13 x 100mm Polypropylene Tube for PE	pkt/1,000
—	4092-330	12 x 75mm Polypropylene Tube	pkt/1,000
B019-3234	4092-400	Centrifuge Tube, 50 ml, Polypropylene, metal free for PE AS-90/91	pkt/500
—	SC0045	4.5ml Polystyrene Tube (10 x 65 mm) for PE AS-90	pkt/500
—	SC0010	Cloverleaf Polypropylene Sample Cup for PE AS-60 and AS-70	pkt/5
—	SC0060	13 x 100mm Polypropylene Tube 6 ml	pkt/1,000
B051-0397	4092-304	1.1ml Cup for PE AS-71/72	pkt/2,000
N9301205		16 x 125mm Polypropylene Tube 16 ml for PE AS-90	pkt/1,000

Note: For higher acidic or organic samples **Polyethylene** is recommended whereas for low acidity aqueous samples **Polystyrene** is suggested

Sundry Items used in Atomic Spectroscopy Instruments ICP, AAS, FIAS, FIMS, MHS etc. and associated Accessories and Auto Samplers etc.

In SEQ. of PERKIN-ELMER Part No.	Part No. CHOICE	Description	Unit of Sale
Quartz Windows			
B006-6549	1000-70	Window, Zeeman 3030, 5X00, HGA-600, 700 Right Hand & HGA-300, 400, 500 Left & Right Hand	each
B011-5576	4097-10	Window, Zeeman 3030, 5X00 & HGA-600, 700 Left Hand	each
B019-8005	4097-13	Window, for THGA Furnace Only	each
Sampling Capillaries and General Use Tubing			
B008-1168	4097-39	Sampling Capillary, AS-1	each
—	4061-01	Teflon Tubing, 3/16 id x 5/16 od	per metre
—	4097-302	Teflon Tubing, .026 id x .050 od	each
—	4097-322	Teflon Tubing, .063 id x .125 od	per metre
B010-1161	4097-35	Sampling Capillary, AS-40	each
B012-9258	4097-37	Sampling Capillary, AS-60, 70, 71, 72	each
0990-8265	4060-38TT	Teflon Tubing, .6 mm id, 10 ft	each



CHOICE Analytical Pty Ltd

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.au

Web : www.choiceanalytical.com.au