

### Peristaltic Pump Tubing:-

- Tygon / PVC
- Coated PVC
- Viton
- Silicone

#### Suitable for use on:-

- Auto-Analyzers
- AAS & ICP Instruments etc.



All tubing features **two colour-coded stops** for easy determination of size and is available in • **Standard Tygon / PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Silicone** • **Santoprene[\*] (Creamy Opaque)**. Tubing overall length is 16", with a 6" centre pumping section separated by two colour-coded stops. Our Viton Tubing (black opaque) has an overall length of 7" and also includes two stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID Inches/mm	Standard Tygon/PVC Pack of 12	Coated PVC Solvent Resistant Tubing Pack of 12	Viton High Chemical & Temperature Resistant Pack of 12	Santoprene Tubing * Pack of 12	Silicone Tubing Pack of 6
Orange/Black	0.005/0.13	AR-4300	AR-4300Y	-	-	-
Orange/Red	0.008/0.19	AR-4301	AR-4301Y	-	-	-
Orange/Blue	0.010/0.25	AR-4302	AR-4302Y	-	AR-4302SA/12	-
Orange/Green	0.015/0.38	AR-4303	AR-4303Y	-	AR-4303SA/12	-
Green/Yellow	0.017/0.44	AR-4304	AR-4304Y	-	-	-
Orange/Yellow	0.020/0.51	AR-4305	AR-4305Y	AR-4305B	AR-4305SA/12	-
White/Yellow	0.022/0.57	AR-4306	AR-4306Y	-	-	-
Orange/White	0.025/0.64	AR-4307	AR-4307Y	AR-4307B	AR-4307SA/12	AR-4307S
Black/Black	0.030/0.76	AR-4308	AR-4308Y	AR-4308B	AR-4308SA/12	AR-4308S
Orange/Orange	0.035/0.89	AR-4309	AR-4309Y	AR-4309B	AR-4309SA/12	AR-4309S
White/Black	0.037/0.95	AR-4310	AR-4310Y	-	-	-
White/White	0.040/1.02	AR-4311	AR-4311Y	AR-4311B	AR-4311SA/12	AR-4311S
White/Red	0.043/1.09	AR-4312	AR-4312Y	-	-	-
Red/Red	0.045/1.14	AR-4313	AR-4313Y	AR-4313B	AR-4313SA/12	AR-4313S
Red/Grey	0.048/1.22	AR-4314	AR-4314Y	-	-	AR-4314S
Grey/Grey	0.051/1.30	AR-4315	AR-4315Y	AR-4315B	AR-4315SA/12	AR-4315S
Yellow/Yellow	0.056/1.42	AR-4316	AR-4316Y	AR-4316B	AR-4316SA/12	AR-4316S
Yellow/Blue	0.060/1.52	AR-4317	AR-4317Y	AR-4317B	AR-4317SA/12	AR-4317S
Blue/Blue	0.065/1.65	AR-4318	AR-4318Y	AR-4318B	AR-4318SA/12	AR-4318S
Blue/Green	0.069/1.75	AR-4319	AR-4319Y	-	-	-
Green/Green	0.073/1.85	AR-4320	AR-4320Y	AR-4320B	AR-4320SA/12	AR-4320S
Purple/Purple	0.081/2.06	AR-4321	AR-4321Y	AR-4321B	AR-4321SA/12	AR-4321S
Purple/Black	0.090/2.29	AR-4322	AR-4322Y	AR-4322B	AR-4322SA/12	AR-4322S
Purple/Orange	0.100/2.54	AR-4323	AR-4323Y	AR-4323B	AR-4323SA/12	AR-4323S
Purple/White	0.110/2.79	AR-4324	AR-4324Y	AR-4324B	AR-4323SA/12	AR-4324S
Black/White	0.125/3.17	AR-4325	AR-4325Y	-	-	-

Many other sizes and types of Tubing and Pump Tubes are also available e.g. **3 stop type**. Please enquire providing full details.

\* Santoprene = PharMed = Ismaprene (All are same material - creamy opaque)



**CHOICE Analytical Pty Ltd**

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : [info@choiceanalytical.com.au](mailto:info@choiceanalytical.com.au)

Web : [www.choiceanalytical.com.au](http://www.choiceanalytical.com.au)