


Listing of Equivalent Items
available from **CHOICE ANALYTICAL** for use on:-



Elementar Vario EL

Consumables Listing for ELEMENTAR VARIO EL

Chemical Reagents and Additives

In sequence of ELEMENTAR Part No.	Part No. CHOICE	Description	Unit of Sale
03 001 119	AEB-4000	Tungstic Oxide, fine powder	10gm
03 679 905	AEB-1002	Copper Oxide, fine wires	50gm
03 679 907	AEB-1008	Quartz (Silica) Chips	50gm
03 679 908	AEB-1017	Quartz (Silica) Wool, very pure, very fine	5gm
03 679 910	AEB-1043	Carbon, high purity, granular	5gm
03 679 914	AEB-1029	Lead Chromate, granular	50gm
03 679 915	AEB-1028	Cerium Oxide, granular	25gm
03 679 924	AEB-1015	Copper Wires, fine, reduced, ampoule	100gm
03 697 900	CA-2174	CA-sorb A (14-20 mesh) 500 gm	500gm
03 960 157	AEB-1044	Silica Gel, granular, indicating	100gm
03 960 352	CA-171	Magnesium Perchlorate, granular	1 lb = 454gm
03 960 355	AEB-1094	Ammonium Chloride	50gm
03 960 356	AEB-1093	Hexamethulene-tetramine	50gm
05 000 086	AEB-1089	Aluminium Oxide Wool	20gm
22 131 365	AEB-1010	Silver Wool, very fine wire	10gm
22 131 365	AEB-1119	Silver Wool, very fine wire unreel	10gm
50 008 467	AEB-1090	Corundum Balls	40gm

Certified Organic Analytical Standards

In sequence of ELEMENTAR Part No.	Part No. CHOICE	Description	Unit of Sale
03 001 118	AEB-2043	Sulphanilic Acid OAS	1gm
03 679 942	AEB-2000	Acetanilide OAS	1gm
03 960 350	AEB-2015	4-Dinitroaniline	1gm
05 000 400	AEB-2004	Benzoic Acid, OAS	1gm

These Organic Analytical Standards are traceable to the National Institute of Standards and Technology.

ctd. page 2



CHOICE Analytical Pty Ltd

ABN 66 052 203 719

Supplies for Laboratory Instruments


Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.au


Web : www.choiceanalytical.com.au

Quartzware & Glassware



In sequence of ELEMENTAR Part No.	Part No. CHOICE	Description	Unit of Sale
05 000 088	AEC-3011	Oxygen Lance, 140 mm	each
05 000 354	AEC-3012	Oxygen Lance, 150 mm	each
1100-1150/4	AEC-3009	Combustion Tube	each
1100-1170/4	AEC-3018	Adsh Finger, silica, 60 mm	each
1100-1171/4	AEC-3016	Support Tube, 133 mm	each
1100-1172/4	AEC-3017	Support Tube, 65 mm	each
1100-1308/4	AEC-3010	Quartz Tube, with ground in ball	each
1100-1317/4	AEC-3021	Ash Finger, silica, 105 mm	each
1100-1331	AEC-3022	Connecting Tube, silica	each
22 133 210/3	AEC-3019	Pyrolysis Tube	each
22 133 211/4	AEC-3020	Pyrolysis Tube Insert, silica	each
22 160 200/4	AEC-3008	Drying Tube	each
22 160 301/4	AEC-3015	U Drying Tube, complete	each

Boats and Capsules




In sequence of ELEMENTAR Part No.	Part No. CHOICE	Description	Unit of Sale
22 137 418	ATD-5022	Boats, tin, 12 x 4 x 4 mm, 21 mg	pk/100
22 137 419	ATD-5023	Boats, tin, 12 x 6 x 6 mm, 38 gm	pk/100
—	ATD-5024	Boats, silver, 12 x 4 x 4 mm	pk/100
22 133 212	ATD-5026	Boats, silver, 12 x 6 x 6 mm	pk/100
03 951 619	ATD-4033	Tin Capsules, smooth, 5.5 x 3.5 mm	pk/100
03 951 620	ATD-4001	Tin Capsules, smooth, 6 x 3 mm	pk/100

cont. page 3



O Rings and Seals



In sequence of

ELEMENTAR
Part No.
Part No.
CHOICE
Description
Unit of Sale

03 002 267	AEE-1075	O Ring, viton	each
03 651 011	AEE-1058	Tubin PTFE 5 x 1 mm x 0.5 mm	each
03 653 103	AEE-1063	O Ring, silicone	each
03 654 603	AEE-1016	O Ring, silicone	pk/10
03 654 627	AEE-1073	O Ring, viton	each
03 654 629	AEE-1066	O Ring, silicone	each
03 654 772	AEE-1017	O Ring, viton	each
03 654 814	AEE-1055	O Ring, silicone	each
03 667 990	AEE-1079	Sealing Ring, PTFE	each
03 667 991	AEE-1080	Sealing Ring, PTFE	each
05 000 043	AEE-1062	O Ring, silicone	each
05 000 092	AEE-1061	O Ring, Kalrez	each
05 000 095	AEE-1074	O Ring, viton	each
05 000 096	AEE-1072	O Ring, viton	each
05 000 140	AEE-1060	O Ring, nitrile	each
05 000 247	AEE-1050	Seal for Oxygen Lance, viton tube, 100 mm	each
05 000 247	AEE-1050	Seal Oxygen Lance, viton, 100 mm	each
05 000 249	AEE-1067	O Ring, viton	each
05 000 353	AEE-1054	Quadring, viton	each
05 000 359	AEE-1070	O Ring, viton	each
05 000 368	AEE-1071	O Ring, viton	each
05 000 369	AEE-1057	O Ring, viton	each
05 000 371	AEE-1076	O Ring, viton	each
05 000 372	AEE-1077	O Ring, viton	each
05 000 382	AEE-1052	O Ring, Kalrez	each
05 000 383	AEE-1051	O Ring, Kalrez	each
05 000 384	AEE-1078	Quadring, viton	each
05 000 385	AEE-1068	O Ring, viton	each
05 000 416	AEE-1069	O Ring, viton	each
1100-1304/4	AEE-2037	Teflon Plug, for reduction tube	each
—	AEC-1020	Cap Kit for threaded scrubbers, caps and teflon washers on silicon seals (AEE-1021)	pk/2

