


Listing of Equivalent Items
available from **CHOICE ANALYTICAL** for use on:-



LECO CHN-2000

Consumables Listing for LECO CHN-2000

LECO Part No.	Part No. CHOICE	Description	Unit of Sale
501-045	AEC-4007	Silica Crucible	each
501-045-125	AEC-4034	Ceramic Crucible	10/pkt
501-059	AR-059	Tin Capsules 14x5mm	100/btl
501-077	CA-077	Iron Chip	908gm=2 1b
501-081	CA-081	Glass Wool	454gm=1lb
501-171	CA-171	Anhydrone	454gm=1lb
501-241	AR-241	Silicone Grease	5.3oz
501-426	AR-426	Comb. Control for Solids	454gm=1lb
501-427	CA-427	Comb. Control for Liquids	454gm=1lb
501-571	AR-571	Copper Capsules 18x5.8mm	100/btl
501-608	AR-1608	Quartz Wool Wadding	200gm
501-609	AR-615	Furnace Reagent	100gm
501-621	AR-621	Copper Turnings	200gm
501-636	AR-636	Vanadium Pentoxide	250gm
501-673	AR-673	Iron Chips	454gm=1lb
502-007	AR-2189	Tube Cleaner 8mmx4ft	each
502-008	AR-570	Tin Plug for AR-571	150/btl
	AR-570T	Plug Insertion Tool	each
502-023	AR-473	Funnel for Reagent Tube	each
502-040	AR-2166	Tin Capsule Large 17x6mm	100/btl
502-040-100	10 x AR-2166	See 502-040 above	
502-042	AR-025	Large Quartz Crucible	each
502-049	CA-049	Reagent N Catalyst	50gm
502-076	AR-571	Copper Capsules	100/btl
502-092	AR-2092	EDTA Standard	50gm
502-141	AEB-4000	Tungsten Oxide Accelerator	10gm
502-167	AR-2167	Tin Capsules XL 20x8mm	100/btl
502-174	AR-168	Alphasolve 20-40 Mesh	250gm
502-177	AR-177	Fine Quartz Wool	28gm=1oz
502-181	AR-2181	Tungstic Oxide, granular	60gm
502-186	CA-2186	Tin Foil Cup	100/btl
502-186-100	10 x CA-2186	see 502-186 above	
502-187	AR-2187	Tin Capsules 9x3.7mm	100/btl
502-188	AR-2188	Aluminium Oxide Pellets	200gm
502-189	AR-189	Copper Sticks Reduced (Ampoule)	100gm
502-190	AEB-1002	Copper Oxide Fine Wires	50gm
502-206	ASD-2003	Silver Capsule 3.5x5mm	100/btl
502-226	AEB-1009	Silver Vanadate	10gm
502-227	ATD-1000	Tin Capsules 4x3.2mm	500/btl

ctd. page 2



CHOICE Analytical Pty Ltd

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.au

Web : www.choiceanalytical.com.au

Consumables Listing for LECO CHN-2000

LECO Part No.	Part No. CHOICE	Description	Unit of Sale
502-231	AR-673	Hi Purity Iron Chip	454gm=1lb
502-295		(smaller pack of 501-621)	
502-304 and -500	AR-2304	Copper Sticks (Jar) <i>We recommend 5xAR-2304 or 5xAR-189 to minimise oxidation after opening and before final use.</i>	100/500gm
502-310	AR-2310	Steel Wool	1lb=454gm
502-311	AR-2311	Furnace Reagent II	100gm
502-338	AR-338	AlphaGel Caps	400/btl
503-032	AR-901	Glass Scoop	each
529-026	AR-026	Crucible (CHN/FP)	50/pkt
529-055	AR-055	Crucible larger size	50/pkt
537-103	AR-103	PolyFlo Tubing	Ft
537-153	AR-153	Ferrules Nylon for AR-103	each
589-551	AR-422	O Ring Piston Seal	each
601-012		O Ring CHN-1000	each
601-133	AR-1133	O Ring	each
601-437	AR-437	Quartz Catalyst Tube	each
601-441	AR-441	Lance Quartz FP-428	each
601-442	AR-442	O Ring Reduction Heater	each
601-482	AR-482	Honeycomb Support	each
601-483	AR-483	Paper (500 sheets)	pack
601-484	AR-484	Ribbon & Spool FP-3/428 FP/CNS-2000 & 932/1000/2000	each
601-490	AR-6490	Foiler Device (use w/502-186)	each
601-504	AR-1504	O Ring Loading Mechanism	each
601-691	AR-1691	CFLEX Tube	1 foot
601-963	AR-1963	Tin Capsule 3.5x4mm Flat Bottom	100/btl
604-541	AEC-4033	Lance	each
605-828	AR-828	O Ring Red/Ox Tube Bottom	each
605-830	AR-830	O Ring Red/Ox Tube Top	each
606-333	AR-6333	O Ring	each
608-343	AR-8343	Heating Element, CHN-1000	each
608-379	CA-8379	Quartz Wool Strips	10/pkt
616-108	AR-6108	Combustion U-Tube CHN-2000	each
616-138	AR-6138	O Ring	each
616-139	AR-6139	O Ring	each
616-146	AR-6146	Lance only	each
616-147	AEE-2199	Lance Assembly Complete	each
760-224	AR-224	O Ring	each
761-484	AR-1000	Printer Ribbon	each
761-747	AR-611	Schutzes Reagent	100gm
761-929	AR-929	Tongs	each
763-265		Smaller pack of 501-081 above (Glass Wool)	
768-980	AR-980	Particle Filter	each
771-030	AR-030	Bubble Filter Tube 4"=102mm	each
772-269	AR-269	O Ring	each
772-520	AR-520	O Ring	each
772-738	AR-738	O Ring	each
772-910	AR-2910	O Ring Piston Seal	each
773-320	AR-320	O Ring	each
773-913	AR-3913	O Ring Loading Head	each
774-425	AR-432	O Ring	each
774-638	AR-638	Disconnect Tube 6"=152mm	each

ctd. page 3


CHOICE Analytical Pty Ltd

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.auWeb : www.choiceanalytical.com.au

Consumables Listing for LECO CHN-2000

LECO Part No.	Part No. CHOICE	Description	Unit of Sale
775-306	AR-5306	Mesh Filter	each
775-314	AR-314	O Ring	each
775-601	AR-601	Reagent Tube (left) 6"=152mm	each
775-880	AR-880	Disconnect Tube 4"=102mm	each
776-279	AR-279	Catalyst Heater Tube	each
776-284	AR-284	O Ring Combustion Tube	each
776-297	AEC-4006	Purification Tube	each
778-114	AR-8114	O Ring	each
778-116	AR-8116	O-Ballast	each
778-120	AR-781	Ballast Gasket	each
780-899	AR-899	Honeycomb Stop	each
783-189	AR-3189	O Ring CHN-1000	each
783-274		Furnace Lance (783-278 = AR-278) plus collar assembly	each
783-278	AR-278	Furnace Lance only, CHN-1000 and CHN-2000	each
783-897	AR-3897	O Ring CHN-1000	each


A complete range of organic standards and tin and silver capsules for microanalysis is available. Please enquire with your requirements

AR-2185

Nitrogen Free Paper Cups

Unit: 100 cups

For use in Nitrogen/Protein analysis. These cups have no Nitrogen and are virtually ashless. They have the same dimensions as the AR-2186 Tin Foil Cups and can be used for dry, solid samples.



AR-1490 (601-490) Foiler Fixture

This easy to use Foiler Fixture is now available through Choice Analytical. This press has a funnel shape design for easy use with the AR-2186 Tin Foil Cups or AR-2185 Nitrogen Free Paper Cups. You may also use any other capsule with a diameter up to 8mm. The sample ends up in a small pellet form, approximately 8mm in diameter.



AlphaGel Caps

The **AR-338** (502-338) AlphaGel Cap is for use on the FP-228, 428 and 528 instruments. These Nitrogen free capsules are for use with dry samples and do not leave a

residue in the crucible. This affordable method of sample containment will save you both time and money. 400 capsules per container.




AR-4961 (614-961) Large Crucible

Large crucible for FP-528 or other applications.

Quick Cap Kit

Choice Analytical is pleased to offer its customers another new product. This helpful Quick Cap Kit **AR-9170** (609-170) will allow you to place dry sample material into the AlphaGel Cap **AR-338** (502-228) for weighing, prior to analysis. The Kit comes complete with one plunger **AR-9171** (609-171) and ten tubes **AR-9172** (609-172). Capsule Holder **AR-4373** (604-373), additional tube and AlphaGel Caps sold separately.



CHOICE Analytical Pty Ltd

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : info@choiceanalytical.com.au

Web : www.choiceanalytical.com.au

ABN 66 052 203 719