

### Peristaltic Pump Tubing:-

- Tygon / PVC
- Coated PVC
- Viton
- Silicone

Suitable for use on:-

**PERKIN-ELMER**

AAS & ICP-OES,  
ICP-MS, FIMS, FIAS Instruments



All tubing features **two or three colour-coded stops** for easy determination of size and is available in • **Standard Tygon / PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Silicone** • **Santoprene (creamy opaque)**. Tubing overall length is 16", with a 6" centre pumping section separated by two colour-coded stops, three stops (72 mm between stops) Our Viton Tubing (black opaque) has an overall length of 7" and also includes two or three stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID inches/mm	Perkin-Elmer Part No.	Standard Tygon/PVC Tubing	Coated PVC Solvent Resistant Tubing	Santoprene Tubing	Pack Size 12/pk	Silicone Tubing	Pack Size 6/pk
<b>TWO-STOP TYPE</b>								
used on ICP-OES systems								
Black/Black	0.030/0.76	0990-8587	116-0549-07	-	-	12/pk	-	-
Black/Black	0.030/0.76	0047-3550	-	116-0533-07	-	12/pk	-	-
Black/Black	0.030/0.76	0047-3552	-	-	-	-	116-0497-08	6/pk
Red/Red	0.045/1.14	N068-0375	116-0549-10	-	-	12/pk	-	-
Red/Red	0.045/1.14	N068-0376	-	116-0533-10	-	12/pk	-	-
Black/White	0.125/3.17	N812-2012	116-0549-21	-	-	12/pk	-	-
Orange/Red	0.008/0.19	N069-5476	116-0549-01	-	-	12/pk	-	-
<b>THREE-STOP TYPE</b>								
used on ICP-MS, FIAS and FIMS systems								
Black/Black/Black	0.030/0.76	B050-6058	SC0056	-	-	12/pk	-	-
White/White/White	0.040/1.02	B050-7692	-	SC0292	-	12/pk	-	-
Red/Red/Red	0.045/1.14	B019-3160	SC0061	-	-	12/pk	-	-
Red/Red/Red	0.045/1.14	B3140730	-	-	SC0310	12/pk	-	-
Yellow/Blue/Yellow	0.060/1.52	B019-3161	SC0065	-	-	12/pk	-	-
Purple/Purple/Purple	0.081/2.06	B019-9034	SC0069	-	-	12/pk	-	-
Purple/White/Purple	0.110/2.79	B3140721	-	-	SC0319	12/pk	-	-
Black/White/Black	0.125/3.17	B050-8310	SC0073	-	-	12/pk	-	-

Many other sizes and types of Tubing and Pump Tubes are also available e.g. Viton. Please enquire with full details of requirement.



**CHOICE Analytical Pty Ltd**

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : [info@choiceanalytical.com.au](mailto:info@choiceanalytical.com.au)

Web : [www.choiceanalytical.com.au](http://www.choiceanalytical.com.au)